Normal Automatic Sequencing

Standard Sequencing



Description

Macrogen Inc. offers DNA sequencing using state-of-the-art robotics and instrumentation to generate high quality sequence data at affordable prices. We have successfully sequenced countless samples of PCR products, plasmids, and BACs from customers all over the world, and our DNA sequencing services are highly rated for their turnaround time, quality of results, and competitive prices.

Features

- 3730xl DNA analyzer
- Real time monitoring from registration to result delivery using LIMS (Laboratory Information Management System)
- Result delivery within 24 hours.
- Performing fundamental analysis service at free of charge (Investigating homology for the protein and DNA sequence level through BlastN or BlastX)

Sample preparation

Template

- Plasmid : Concentration 100ng/µl, Volume 20µl
- PCR product (purified) : Concentration 50ng/µl, Volume 20µl
- PCR product (non-purified) : Concentration 100ng/μl, Volume 20~30μl
 Example) In case of 300bp of PCR product, if you send 30ng/μl of a sample then better results are obtained(100bp≒10ng/μl). It is recommended to check the PCR product concentration through electrophoresis rather than Nano-drop.

Primer

- Concentration: 5 pmole/µl (5pM), Volume: 20µl
 Example) For 5 Reactions: 20µl, as the reaction number increases increment 2µl each.
- Please separate template and primer when you send primers.
- Primer Tm : 55~60°C, GC content : 40%~60%.

Order & Result

 Based on normal plasmid template basis, guarantee phred score more than 20, read-length more than 650 bp



*Plate Sequencing

Description

For the customers with large number of samples processed in 96-well plates, this will be a better option. "Plate Sequencing" refers to the orders charged per one entire 96-well plate. Samples can be sent as clones, un-purified PCR products and purified DNA.

Benefits

- Low Cost & Best Price
- High Quality Service
- Quick Turnaround
- Shipping Free

5 6

 A
 1
 9
 17
 25
 33
 41
 49
 57
 65
 73
 81
 89

 B
 2
 10
 18
 26
 34
 42
 50
 58
 66
 74
 82
 90

C 3 11 19 27 35 43 51 59 67 75 83 91

 D
 4
 12
 20
 28
 36
 44
 52
 60
 68
 76
 84
 92

 E
 5
 13
 21
 29
 37
 45
 53
 61
 69
 77
 85
 93

 F
 6
 14
 22
 30
 38
 46
 54
 62
 70
 78
 86
 94

 G
 7
 15
 23
 31
 39
 47
 55
 63
 71
 79
 87
 95

H 8 16 24 32 40 48 56 64 72 80 88 96

7 8 9 10 11

Sample preparation

Easy order system

- 96 Well Plate:
- Please number your samples as shown the well box (top to bottom, left to right, A1 to H12):

Service Type



- One single primer should be applied to the whole plate, so that we can simply put the primers to the corresponding sample plate.
- We will provide the results as the following file name.

[Plate_name]-[Well_position]_[Primer_name].ab1

2primer / 1plate

 Only two custom primers (or two universal primers) in tubes could be acceptable for applying this Plate sequencing.

*Note : For three or more primers, your order will be regarded as Single sequencing reactions

Multiple Primers / 1plate

 "Multiple Primers / 1 Plate can be applicable in case you send primers in the plate which has same well position information with the sample plate regardless of primer's quantity.
 *Note : Each paired primer should be filled in the 96-well plate which has same well position with the paired sample plate

Value Service

- Data registration and the process can be monitored in real time. For customer's convenience, the result can be delivered as downloadable file (ab1, txt, pdf, phd format) through your LIMS and e-mail.
- Also, Macrogen lifetime storage service will be stored all of your result files (ab1 format) free of charge up to fifteen (15) years from the date of data creation.
- All process from registration to result sending can be done via on-line sequencing system (http://dna.macrogen.com).



MACROŬGEN

 Macrogen Inc.
 #60-24, Gasan-dong, Geumchen-gu, Seoul 158-781, Korea

 Macrogen Europe
 Meibergdreef 39, 1105 AZ, Amsterdam, the Netherland

 Macrogen USA
 9700 Great Seneca Hwy, Rockville, MD 20850, United States

 Macrogen Japan
 1-1-1 Sakuragaoka, Setagaya-ku, Tokyo, 156-8502, Japan

e-mail : info@macrogen.com e-mail : europe@macrogen.com e-mail : customerservice@macrogenusa.com e-mail : customer@macrogen-japan.co.jp